

## Letter: Beyond frame expansion: interpreting the implications of routine post-dilatation

Matteo Sturla, MD; Pier Pasquale Leone, MD; Azeem Latib\*, MD

\*Corresponding author: Interventional Cardiology, Montefiore-Einstein Center for Heart and Vascular Care, Montefiore Medical Center/Albert Einstein College of Medicine, 111 E 210<sup>th</sup> St, Bronx, NY, 10467, USA. E-mail: alatib@gmail.com

We read with great interest the study by Husain et al<sup>1</sup> and applaud the authors for this elegant and methodologically rigorous investigation. The authors combined benchtop and *in vivo* evaluations of routine “double-tap” post-dilatation at nominal volume following transcatheter aortic valve implantation with a balloon-expandable valve (BEV). Using a fluoroscopic calibration of valve commissural post height, they demonstrated that nominal-pressure reinflation increased frame expansion by approximately 9% across all SAPIEN 3 Ultra (Edwards Lifesciences) valve sizes, without excess procedural or 30-day adverse events<sup>1</sup>.

In a study from our group (Leone et al<sup>2</sup>), we evaluated a similar concept solely focusing on the 23 mm SAPIEN 3 Ultra, a size particularly susceptible to clinically relevant underexpansion and prosthesis-patient mismatch (PPM). Quantitative frame assessment was performed using pigtail-based fluoroscopic calibration, referencing the 10 mm spacing between radiopaque markers. This method provides a wider calibration baseline and minimises proportional measurement error compared with catheter width-based methods<sup>3</sup>. Nonetheless, this technique is inherently subject to interobserver variability, and unlike the commissural post height method, no correlation with computed tomography is available<sup>4</sup>. Distinct from the DOUBLE-TAP protocol, our trial protocol allowed additional volume of the delivery balloon during the “double-tap” (+1-2 mL in 37% of patients) and selective use of non-compliant balloons. Compared with a control group without routine post-dilatation, this strategy yielded greater valve expansion, lower transvalvular gradients, and a reduced incidence of PPM. The haemodynamic benefit

was primarily observed in patients with normal flow, while those with low-flow physiology derived little measurable improvement.

Both studies, albeit with different methodologies, converge on the principle that optimisation of frame geometry may favourably influence forward-flow haemodynamics. This effect is likely most pronounced in smaller valves and in patients with a larger body surface area, where even modest gains in frame expansion can improve transvalvular gradients and attenuate PPM (non-structural valve dysfunction)<sup>5</sup>. By contrast, in larger valve sizes, post-dilatation may have little measurable impact on immediate haemodynamics yet could still mitigate leaflet malcoaptation and pinwheeling related to underexpansion – phenomena that can exacerbate leaflet stress, induce asymmetric coaptation, and predispose patients to hypoattenuated leaflet thickening over time; these mechanisms are plausibly linked to structural valve degeneration and impaired durability<sup>6,7</sup>. Taken together, these considerations reinforce the idea that not all valves or patient anatomies require routine post-dilatation for immediate haemodynamic optimisation, while the potential long-term benefit on leaflet durability remains an open and important question for future study.

Nevertheless, direct evidence linking optimised expansion to improved durability remains inferential. The recent PARTNER 3 data demonstrated that, despite post-dilatation in only ~21% of cases, the mean gradient at 7 years remained stable at 12 mmHg, and the cumulative reintervention rate was 6% between years 1 and 7<sup>8</sup>. Although these data imply that routine post-dilatation may not be essential for maintaining favourable long-term valve function and durability, both

our data and the DOUBLE-TAP findings confirm that there is little to no price to pay in terms of adverse events, thus, a tailored yet liberal approach could be reasonable to adopt.

While we commend the authors for highlighting the prevalence of BEV underexpansion, the clinical significance of nominal-volume “double-tapping” remains to be proven beyond geometric correction. Extended durability will be essential to determine whether routine post-dilatation should evolve from an optimisation tool into a standardised procedural step.

## Authors' affiliation

Montefiore-Einstein Center for Heart and Vascular Care, Montefiore Medical Center, Albert Einstein College of Medicine, Bronx, NY, USA

## Conflict of interest statement

A. Latib has served on the advisory board for Medtronic, Abbott, Boston Scientific, Edwards Lifesciences, Shifamed, NeoChord, V-dyne, and Philips. The other authors have no conflicts of interest to declare relevant to the contents of this paper.

## References

1. Husain A, Jelisejevas J, Millar K, Alnuwakhtha A, Ferkh A, Delarive J, Gill H, Altheekrallah J, Offen S, Tzimas G, Ye J, Leipsic JA, Blanke P, Sathananthan J, Wood DA, Sellers SL, Meier D, Webb JG. Routine post-dilatation at nominal volume to optimise the expansion of balloon-expandable valves: the DOUBLE-TAP study. *EuroIntervention*. 2025; 21:e1159-68.
2. Leone PP, Sturla M, Spring AM, Echarte-Morales J, Scotti A, Ludwig S, Coisne A, Slipczuk L, Assafin M, Chau M, Ho EC, Granada JF, Latib A. Routine postdilation after 23 mm Sapien 3 Ultra implantation in the aortic position. *Catheter Cardiovasc Interv*. 2024;104:559-69.
3. Jacobs JH Jr, Bove AA, Smith HC, Chesebro JH. Use of a metal ring-marked catheter for geometric calibration in quantitative coronary angiography. *Cathet Cardiovasc Diagn*. 1988;15:121-4.
4. Husain A, Tzimas G, Akodad M, Jelisejevas J, Offen S, Millar K, Leipsic JA, Blanke P, Wood DA, Sellers SL, Webb JG, Meier D, Sathananthan J. Novel Real-Time Fluoroscopic Assessment Method of THV Expansion Following Balloon-Expandable TAVR. *Circ Cardiovasc Interv*. 2025;18:e014617.
5. Leone PP, Regazzoli D, Pagnesi M, Cannata F, Mangieri A, Hokken TW, Costa G, Barbanti M, Teles RC, Adamo M, Taramasso M, Reifart J, De Marco F, Giannini F, Kargoli F, Ohno Y, Saia F, Buono A, Ielasi A, Pighi M, Chiarito M, Bongiovanni D, Cozzi O, Stefanini G, Ribichini FL, Maffeo D, Chizzola G, Bedogni F, Kim WK, Maisano F, Tamburino C, Van Mieghem NM, Colombo A, Reimers B, Latib A. Implantation of contemporary transcatheter aortic valves in small aortic annuli: the international multicentre TAVI-SMALL 2 registry. *EuroIntervention*. 2023;19:256-66.
6. Fukui M, Bapat VN, Garcia S, Dworak MW, Hashimoto G, Sato H, Gössl M, Enriquez-Sarano M, Lesser JR, Cavalcante JL, Sorajja P. Deformation of Transcatheter Aortic Valve Prostheses: Implications for Hypoattenuating Leaflet Thickening and Clinical Outcomes. *Circulation*. 2022;146:480-93.
7. Fuchs A, De Backer O, Brooks M, de Knecht MC, Bieliuskas G, Yamamoto M, Yanagisawa R, Hayashida K, Søndergaard L, Kofoed KF. Subclinical leaflet thickening and stent frame geometry in self-expanding transcatheter heart valves. *EuroIntervention*. 2017;13:e1067-75.
8. Leon MB, Mack MJ, Pibarot P, Hahn RT, Thourani VH, Kodali SH, Généreux P, Kapadia SR, Cohen DJ, Pocock SJ, Zhang Y, Szerlip M, Ternacle J, Malaisrie SC, Herrmann HC, Szeto WY, Russo MJ, Babaliaros V, Nazif T, Webb JG, Makkar RR; PARTNER 3 Investigators. Transcatheter or Surgical Aortic-Valve Replacement in Low-Risk Patients at 7 Years. *N Engl J Med*. 2025 Oct 27. [Epub ahead of print].